

# Pull-Tab Tin Can Opener

## USER QUICK GUIDE

### Introduction

The Pull-Tab Tin Can Opener is a device to help those with arthritis, limited finger dexterity or limited finger strength open tin cans with pull-tabs.

### Features



**Pull-Tab Hook:** The point which grabs the Pull-Tab so it can be opened.

**Handle:** Where the User holds the device.

### Usage

#### Step 1

Slide the hook under the pull-tab of the can.



#### Step 2

Slide the hook into the pull-tab so that the outer ring sits into the groove and the seal of the can breaks.

# Pull-Tab Tin Can Opener

## USER QUICK GUIDE



### Step 3

Lever the can opener backwards to pull open the top of the can.



### Step 4

Continue to pull the can opener and can lid off until the can is open.

# Pull-Tab Tin Can Opener

## USER QUICK GUIDE



### Compatibility

This device is compatible with any tin can with a pull-tab, such as tuna, pet food or soup cans.

### Specifications

This device is 17.7 cm long, 7.5 cm wide, and 1 cm tall.

### Cleaning

Use sterilizing wipe or spray. Do not place in dishwasher or expose to high heat.