

Pull-Tab Tin Can Opener

USER QUICK GUIDE

Introduction

The Pull-Tab Tin Can Opener is a device to help those with arthritis, limited finger dexterity or limited finger strength open tin cans with pull-tabs.

Features



Usage

Step 1

Slide the hook under the pull-tab of the can.



Step 2

Slide the hook into the pull-tab so that the outer ring sits into the groove and the seal of the can breaks.



© 2022 by Neil Squire / Makers Making Change.

This work is licensed under the CC BY SA 4.0 License: <http://creativecommons.org/licenses/by-sa/4.0>

Files available at <https://makersmakingchange.com/project/pull-tab-tin-can-opener/>

Pull-Tab Tin Can Opener

USER QUICK GUIDE

Step 3

Lever the can opener backwards to pull open the top of the can.



Step 4

Continue to pull the can opener and can lid off until the can is open.



Compatibility

This device is compatible with any tin can with a pull-tab, such as tuna, pet food or soup cans.

Specifications

This device is 17.7 cm long, 7.5 cm wide, and 1 cm tall.

Cleaning

Use sterilizing wipe or spray. Do not place in dishwasher or expose to high heat.